

## OVERSIZE BODY BOX ASSEMBLY INSTRUCTIONS

(No tools or adhesives required)

The container is designed to simplify the cremation process. The container is cut and folded so that one person may load and close it. Each unit includes components for two cremations.

Contents are: 2 covers, 2 containers, 2 plastic body bags, 4 straps, 2 liners, and 4 clips. The process starts with the container opened up and flat on a table and goes as follows:



1

Photo of box.



2

Box with lid open.



3

Body box containers.



4

Cut straps on container box.




5

Container ready to fold.




6

The plastic bag is for encapsulating the body so the fluids don't leak. The bands are for closing the box.



7

Photo of every part labeled.




8

Grab the square end panel and bring panels B & C upright. Fold panel C into width of box.


*“Leading the Change in the Funeral Industry”*

## OVERSIZE BODY BOX (CONTINUED...)




9

Grab panel E, bring upright, placing panel D against C.




10

Fold panel E over and against C, folding panel F inside panel B.



11

Handholes on panel F and B will line up. Push tab into F to lock in place. (repeat this process on the opposite end and then load the body)



12

Fold both panels H so that they are upright and then fold up panel I, placing both panels H between E and D. Push tab into G to lock in place.



13

To insert the liner in the box, fold outside panels so they are upright, then insert the liner inside the container pushing down to the bottom. (This will achieve the weight Capacity)



14

Cover container with lid. Use enclosed straps if necessary.

*“Leading the Change in the Funeral Industry”*